



OXYGEN CONCENTRATOR

- France imported CECA molecular sieve
- High oxygen purity $93\% \pm 3\%$
- Long life span more than 10000 hours
- Suitable for 24 hours continuously work
- 2 years warranty
- CE ISO approved
- OEM available

 **live** / CARE FOR YOU!
www.zzolive.com

SPECIFICATION

Model	OLV-3	OLV-3W	OLV-5	OLV-5W
Purity Oxygen	≥93%±3 (V/V)			
Voltage	AC220V±22V, 50Hz±1Hz/AC110V±15V, 60Hz±1Hz			
Outlet Pressure	20kPa~60kPa			
Flow Rate	0-3L/min		0-5L/min	
Sound Level	less than 45db			
Atomization		≥0.15mL/min		≥0.15mL/min
Power Consumption	260Watts		300Watts	
Net Weight	19kg		22kg	
Dimension	490×270×560 (mm)			
Warranty	2 years			
Standard Accessories	1× Power cord 1× Mouth-contained nebulizer (Only atomization models) 1× Spare connecting tube of humidifier bottle 1× Spare dust filter 3× Spare inner filter 2× Nasal pipe 1× user manual 1× warranty card			

CERTIFICATE

